High voltage power supply

PF1055

Quad output, High Voltage Rack Mount Power Supply designed for use with PHOTONIS bipolar time-of-flight detectors



Applications

- ✓ Channeltron®
- ✓ Electron multipliers
- Microchannel plate detectors
- Discrete dynode detectors
- ✓ Photomultiplier tubes
- Mass spectrometry
- Analytical instrumentation
- Geiger tubes
- Proportional counters

Features

- Input and output, scintillator and PMT supply provides drive voltages for MCP
- ✓ Four panel mounted output voltage display
- ✓ Individual front panel voltage adjustments
- Supply adjustment for positive and negative polarity
- √ 90-240 V ac input voltage
- Overload and short-circuit protected
- ✓ SHV output connectors
- ✓ All solid-state design
- ✓ Exceptional stability
- ✓ Low ripple < 0.05%</p>
- Compact design
- Very economical

Model PF1055 is a solid-state, rack mount, high voltage power supply providing the four DC output voltages required to operate PHOTONIS Bipolar Time-Of-Flight detectors. The polarities can be set to tune the detector to detect either positive or negative ions. This unique power supply is specifically designed to provide over-voltage and arc protection for the sensitive components of the Bipolar TOF detectors. This product eliminates the supply synchronization problems associated with the use of individual power supplies. The PF1055 offers high precision and stability in a compact, economical package which makes it ideal for many laboratory applications.

Model PF1055 Quad Output Power Supply

Specifications	
Dimensions	17" (D) x 19" (W) x 5.25" (H)
Input Voltage	90 to 240 V ac, 50-60 Hz
PMT supply (referenced to ground)	0 to 1 kV DC, 2 mA
Scintillator supply (referenced to MCP Supply)	0 to 3 kV, 100 μA
Bias Voltage (referenced to ground)	0 to 10 kV DC, 200 μA
MCP Supply (referenced to Bias Voltage)	0 to 1 kV DC, 100 μA
Output Polarity (Scintillator, Bias, MCP)	Positive or Negative
Ripple	< 0.05%
Stability	< 0.05% / 8 hours (PMT Supply, Bias Voltage)
Line Regulation	0.05%
Load Regulation	0.05%
Output connectors	SHV male (PMT Kings type 1704-1 or equiv., Scintillator, Bias, MCP Kings type 1064-1 or equiv.)



The PF1055 is specifically designed to operate with the PHOTONIS Bipolar Time-of-Flight detectors



The information furnished is believed to be accurate and reliable, but is not guaranteed and is subject to change without notice. No liability is assumed by PHOTONIS for its use. Performance data represents typical characteristics as individual product performance may vary. Customers should verify that they have the most current PHOTONIS product information before placing orders. No claims or warranties are made as to the application of PHOTONIS products. Pictures may not be considered as contractually binding. This document may not be reproduced, in whole or in part, without the prior written consent of PHOTONIS.

04/11/2012