

Hyper-Cam Mini



*IR HYPERSPSCTRAL
IMAGING SYSTEM*

KEY FEATURES



**LOW SWAP, MODULAR DESIGN
INCLUDING OPTICAL HEAD (OH) AND
CONTROL AND PROCESSING BOX (CPB)**



**FOURIER-TRANSFORM IMAGING
SPECTROSCOPY CAPABILITY WITH
USER-SELECTABLE SPECTRAL
RESOLUTION UP TO 4 cm^{-1}**



**320 X 256 PIXEL COOLED SLS DETECTOR
SENSITIVE OVER MWIR ($2.9 - 5.2\text{ }\mu\text{m}$) OR
LWIR ($7.5 - 12.4\text{ }\mu\text{m}$) SPECTRAL RANGE**



**CALIBRATED SPECTRAL RADIANCE
DATA PRODUCT**

The Hyper-Cam Mini is an advanced, compact infrared hyperspectral imaging system that combines high spatial, spectral, and temporal resolution capabilities. Sensitive in the MWIR ($2.9 - 5.2\text{ }\mu\text{m}$) or LWIR ($7.5 - 12.4\text{ }\mu\text{m}$) spectral range, the Hyper-Cam Mini is well suited for the analysis of a broad range of gas, mineral, and other target materials. Reduced size, weight, and power specifications ensure that the Hyper-Cam Mini can be deployed into even the most difficult-to-access field locations by a single operator.

Hyper-Cam Mini



Industrial gas detection & identification



Reveal Pro 6 full-featured scientific software

SPECIFICATIONS	Hyper-Cam Mini MWF	Hyper-Cam Mini xLW
Detector Type	Cooled SLS	Cooled SLS
Detector Format	320 x 256 pixels	320 x 256 pixels
Spectral Range	2.9 – 5.2 μm (1920 – 3450 cm^{-1})	7.5 – 12.4 μm (806 – 1333 cm^{-1})
Field of View	14° x 11°	14° x 11°
Maximum spectral resolution	0.5 cm^{-1}	0.5 cm^{-1}
Noise Equivalent Spectral Radiance (typical)	< 10 $\text{nW/cm}^2.\text{sr.cm}^{-1}$	< 30 $\text{nW/cm}^2.\text{sr.cm}^{-1}$
Radiometric Accuracy	< 2 K	< 3 K
Dimensions	20 x 27 x 21 cm (OH), 21 x 21 x 22 cm (CPB)	20 x 27 x 21 cm (OH), 21 x 21 x 22 cm (CPB)
Weight	< 8.2 kg (OH), < 4.2 kg (CPB)	< 8.2 kg (OH), < 4.2 kg (CPB)
Power Consumption	< 410 W (max), < 270 W (steady-state)	< 410 W (max), < 270 W (steady-state)
Operational Temperature	-10 °C to +50 °C	-10 °C to +50 °C
Storage Temperature	-20 °C to +60 °C	-20 °C to +60 °C

sales@telops.com



exosens.com

EXOSENS
REVEAL THE INVISIBLE

© Telops. The information furnished is believed to be accurate and reliable, but is not guaranteed and is subject to change without notice. No liability is assumed by Telops group of companies nor by any Exosens Group companies. Performance data represents typical characteristics as individual product performance may vary. Customers should verify that they have the most current product information from the Telops group of companies before placing orders. Texts and pictures may not be considered as contractually binding. This document may not be reproduced, in whole or in part, without the prior written consent of Telops.